

carumed

Because we care!

Otoscreen 2 video otoscope for ear canal and drum



The Otoscreen 2 video otoscope is a portable, user-friendly device used to examine the ear canal and ear drum. It can be used in hospitals or medical practices.

The Otoscreen 2 video otoscope means you can perform examinations more comfortably and refine your diagnosis by giving you the perspective necessary for a detailed examination.

carumed

Företagshusvägen 5, 244 93 Kävlinge, Sweden

Phone: +46 (0) 46-73 41 00

info@carumed.se, www.carumed.se

carumed

Because we care!

Otoscreen 2 Specifications



Display

- screen colour, high resolution LCD
- deimensions diagonally 60 mm (2.4 inches)
- resolution 320 x 240 pixels
- charge indicator light low battery

Materials

- camera head stainless steel 18/8
- body ABS
- screen protector polycarbonate

Lighting

- technology LED
- lamps 4

Controls

- button on/off + image capture function
- socket RCA video (NTSC), cable supplied

Power

- source 1.5V LR6 AA batteries (supplied)
- intensity 10mA

Protection

- against electric shock type BF
- against water splashes IPX0

Dimensions and weight

- dimensions 75 (L) × 63 (W) × 212 (H) mm
- weight approx. approx. 180 g (without the batteries)

Temperature

- during operation +10 à +40° C (30 à 75 %, without condensation)
- during storage -20 à +49° C (≤ 95 %, without condensation)

Atmospheric pressure

- during operation 70 kPa à 106 kPa

- Speculum rigid
- use single-use
- material PP (PolyproPylene)

Applicable standards

- EN/IEC 60601-1 ; EN/IEC 60601-1

Warranty

- 1 YEAR
- Excluding wear parts, batteries and speculums.
- (under normal usage conditions)
- Legal warranty LIFETIME
- against all latent manufacturing faults and defects
- from the date of purchase.

Accessories

- Single-use rigid ear speculums, 5 mm (box of 250)
- 7,5 cm long camera
- 15 cm long camera
- 30 cm long camera

carumed

Företagshusvägen 5, 244 93 Kävlinge, Sweden

Phone: +46 (0) 46-73 41 00

info@carumed.se, www.carumed.se