

# CARUMED

*Because we care!*

## Basic Plus Disposable Latex Gloves



### Disposable Latex Gloves

- Unsterile
- AQL: 1,5
- Fulfills EN455 and EN420
- CE mark
- White
- 100pcs per box
- 10 boxes per carton
- 240mm long

#### Art.no Prepowdered

201000	Small
201001	Medium
201002	Large
201003	Extra Large

#### Art.no Powderfree

201004	Extra-Small
201005	Small
201006	Medium
201007	Large
201008	Extra Large

**CARUMED**

Företagshusvägen 5, 244 93 Kävlinge, Sweden

Tel: +46 (0) 46 - 73 41 00

info@carumed.se, www.carumed.se