

CARUMED

Because we care!

Blue Basic Plus Powderfree Nitrile Gloves



Powderfree Nitrile Glove with a good tactile sensitivity and high wearing comfort.
Very elastic and versatile.

- Nitrile
- Powderfree
- Non-sterile
- Grippy
- AQL 1,5
- CE mark
- Fulfills EN 455, EN 420 and EN 374
- Length: approx 240mm
- Blue
- 200pcs/box, 10boxes/ carton

Art.no		202031	Medium
202029	Extra-Small	202032	Large
202030	Small	202033	Extra-Large

CARUMED

Företagshusvägen 5, 244 93 Kävlinge, Sweden

Tel: +46 (0) 46 - 73 41 00

info@carumed.se, www.carumed.se