

# CARUMED

*Because we care!*

## Mint Basic Plus

### Powderfree Nitrile Gloves



Powderfree Nitrile Glove with a good tactile sensitivity and high wearing comfort.  
Very elastic and versatile.

- Nitrile
- Powderfree
- Non-sterile
- Grippy
- AQL 1,5
- CE mark
- Fulfills EN 455, EN 420 and EN 374
- Length approx. 240mm
- Mint color
- 200pcs/box, 10boxes/ carton

**Art.no**

202034 Small  
202035 Medium  
202036 Large  
202037 Extra-Large

---

**CARUMED**

Företagshusvägen 5, 244 93 Kävlinge, Sweden

Tel: +46 (0) 46 - 73 41 00

info@carumed.se, www.carumed.se