

carumed

Because we care!

Polymer Plus

Polymer coated Latex Gloves



High quality latex gloves, powder free, polymer coated inside and chlorinated outside.

- Latex
- Powder free
- AQL 1,5
- CE mark
- Bright beige
- Fulfills EN455 and EN420
- 100 pcs per box
- 10 boxes per carton

Art.no 201015 Extra-Small
201016 Small
201017 Medium
201018 Large
201019 Extra-Large

carumed

Företagshusvägen 5, 244 93 Kävlinge, Sweden

Tel: +46 (0) 46 - 73 41 00

info@carumed.se, www.carumed.se