

# carumed

*Because we care!*

## Style Black Powderfree Nitrile Gloves



Black, powderfree nitrile gloves with good grip. These gloves are specially popular with hairdressers, tattoo artists and some restaurants due to its black color. This glove has everything you demand from a nitrile glove at the same time it has an appealing color.

- Nitrile
- Powderfree
- Non-sterile
- Grippy
- AQL 1,5
- CE mark
- Fulfills EN 455 and EN 420
- Length: approx. 240mm
- Black
- 100pcs/box, 10boxes/  
carton

#### Art.no

- 202006 Small
- 202007 Medium
- 202008 Large
- 202009 Extra-Large
- 202010 Extra-Extra-Large

---

## carumed

Företagshusvägen 5, 244 93 Kävlinge, Sweden

Tel: +46 (0) 46 - 73 41 00

info@carumed.se, www.carumed.se