



Tablett box with alarm

This tablett box with alarm and a push button opening helps you keep your tablet schedule throughout the day. Tablett box has a simple, intuitive programming, alarm functionality for missed tablets, convenient backlight that makes it a great help when planning today's dietary supplements or medications.



You can program the Electronic Tablett Box with up to 7 times per day. When the alarm goes off, the display will flash, and you can program it to either alert or vibrate or both.

Alarm timer on tablettdispensern can be removed and worn separately, if desired. The unit has a secure locking mechanism but is easily opened by pressing a button.

Dimensions: 15,2 x 10,2 x 3,1 cm
Color: White/blue

Article number: 120240



Företagshusvägen 5, 244 93 Kävlinge, Sweden
Phone: +46 (0) 46-73 41 00
info@carumed.se, www.carumed.se