

CARUMED

Because we care!

Module Tray of Plexi Glass

Small

Art.no. 110235



Tray of plexi glass (acrylic) with lose dividers (included). Easy to adjust for different needs. Can be used separately or together with for instance with the larger module tray (110230). The size is suitable for preparing for instance needles, band aids or similar. For use outside the care field, it can also be used to organize smaller belongings. Easy to clean.

Length * Width: 22,5*8,5cm
Height edge: 3,5cm
Height handle: 5,5cm

CARUMED

Företagshusvägen 5, 244 93 Kävlinge, Sweden

Tel: +46 (0) 46 - 73 41 00

info@carumed.se, www.carumed.se