

carumed

Because we care!

Sharp Disposable Protective Lid of polycarbonate

Lid by 3mm uncrushable polycarbonate suited for cardboard sharp disposable containers with an opening of 40*40cm (55 and 110 liters). Prevents contagious substances to come out in the room.

The lid is transparent which allows easy visual access to see when container change is due.



Art.No. 172002

carumed

Företagshusvägen 5, 244 93 Kävlinge, Sweden

Tel: +46 (0) 46 - 73 41 00

info@carumed.se, www.carumed.se