

CARUMED

Because we care!

Infrared Fever Thermometer SC 53 FH



The thermometer SC 53 FH can be used to measure body temperature in the ear or over the forehead. The thermometer has 10 memory spaces for measured values.

The measuring only takes two seconds and a successful reading is advised by a beep. SC 53 FH is delivered with a table stand. The unit can be placed in this stand when it's not being used. When it is in this stand by mode, the time and the room temperature will be shown on the display.

Measuring in the ear is only recommended for persons older than 6 months (as the ear canals are not fully developed until then). Always measure in the same ear to make sure the same conditions apply.

Measuring over the forehead is conducted by moving the thermometer from one end of the forehead to the other.

Art.No. 240320



Type	4-i-1: ear thermometer, forehead thermometer room temperature thermometer and clock
Measure time, ear	2 sec
Measure time, forehead	2 sec
Water tight	No
Measure range	10,0°C - 50°C
Resolution	0,1°C
Precision	±0,2°C at 36,1°C-39°C ±0,3°C at 10,0°C-36,0°C ±0,3°C at 39,1°C till 50,0°C
Memory	10 measure values
Storage conditions	10-50°C, humidity <95%
Working conditions	18°C-28°C, humidity <95%
Battery indicator	Replace battery when "empty battery" symbol appears
Battery	1 * CR 2032, 3V
Dimensions	110 * 33 * 22 mm
Weight	ca 42 g (with battery)

CARUMED

Företagshusvägen 5, 244 93 Kävlinge, Sweden

Tel: +46 (0) 46 - 73 41 00

info@carumed.se, www.carumed.se