



Because we care!

Ear Speculum



Ear Speculum, 2,5 mm
Art.No. 240200



Ear Speculum, 4 mm
Art.No. 240201

Disposable ear speculums, made by recycled PP (free from BPA, phthalates, latex and PVC). Available in sizes 2,5 and 4 mm. Disposable speculums are hygienic as they prevent cross contamination. Secure fit, no sharp edges. Black. CE marked as class I product.

250 per box, 30 boxes per carton (7.500 pcs).

Fits the most common otoscopes on the market as for example:

Heine Alpha+, Heine K180, Heine BETA 400, 800, Heine MINI 2000 + FO, Heine MINI 3000+FO+3000, Heine BETA 200, KaWe Piccolight C, KaWe Combilight F.O., KaWe Eurolight F.O., Riester Ri-mini, Riester Ri Scope L1, Riester Ri Scope L2 and Riester E-scope.



Företagshusvägen 5, 244 93 Kävlinge, Sweden
Tel: +46 (0) 46 - 73 41 00
info@carumed.se, www.carumed.se