



Because we care!

Ear Speculum Eco+



Ear Speculum Eco+, 2,5 mm
Art.No. 240203



Ear Speculum Eco+, 4 mm
Art.No. 240204

Disposable ear speculums, PP made from renewable raw material (free from BPA, phthalates, latex and PVC). Available in sizes 2,5 and 4 mm. Disposable speculums are hygienic as they prevent cross contamination. Secure fit, no sharp edges. Black. CE marked as class I product. The production of the Eco+ speculums is certified by the independant organisation ISCC and annual audits are being performed.

250 per box, 30 boxes per carton (7.500 pcs).

Fits the most common otoscopes on the market as for example:
Heine Alpha+, Heine K180, Heine BETA 400, 800, Heine MINI 2000 + FO, Heine MINI 3000+FO+3000, Heine BETA 200, KaWe Piccolight C, KaWe Combilight F.O., KaWe Eurolight F.O., Riester Ri-mini, Riester Ri Scope L1, Riester Ri Scope L2 and Riester E-scope.



Företagshusvägen 5, 244 93 Kävlinge, Sweden
Tel: +46 (0) 46 - 73 41 00
info@carumed.se, www.carumed.se