

carumed

Because we care!

Digital Fever thermometer DT-12B



- Waterproof
- Mercury free
- With acoustic signal
- Auto switch-off after 10 minuter
- Measuring range: 32-42°C
- Accuracy: $\pm 0,05^{\circ}\text{C}$
- LCD-display
- LR-41 (button cell battery 1,5V)
- Dimensions: 129 * 22 * 14 mm
- Delivered in transparent box
- Medical product according to class Im

Art.No. 240314

carumed

Företagshusvägen 5, 244 93 Kävlinge, Sweden

Tel: +46 (0) 46 - 73 41 00

info@carumed.se, www.carumed.se