

# CARUMED

*Because we care!*

## Infrared Fever Thermometer SC 8280

For professional use





Professional infrared thermometer. To be used 3-8 cm from the forehead or object that shall be measured. Measuring time approximately 1 second. Possible to use for measuring the temperature of objects between 0 to 90°C.

Mainly designed for hospital and professional use, but also possible to use in home environment.

One thermometer per box.

**Art.No. 240330**



CE 0197

Type	Professional infrared thermometer
Measuring time	~1 sec
Measure interval, forehead	32,2°C to 43,3°C
Measure interval, objects	0,0°C to 90,0°C
Unit	°C
Resolution	0,1°C
Precision	±0,3°C at 32,2°C to 43,0°C
Operational temperature	±0,4°C at all other temperatures within the interval 0,0°C to 90,0°C
Memory	+16°C to +35°C at 15 to 90% rel. humidity (non condensed form)
Lifetime batteries	10 measure values
Battery	~5.000 measures with LR03 batteries
Auto shutoff	2 * AAA (R03 or LR03)
Weight	After ~60 sec
Certifications and standards	ca 132 g (with battery)
	ISO 13485, ISO 80601-2-56, CE 0197

## CARUMED

Företagshusvägen 5, 244 93 Kävlinge, Sweden

Tel: +46 (0) 46 - 73 41 00

info@carumed.se, www.carumed.se