

CARUMED

Because we care!

Schine Pill Box



Small, red
Art.No. 120045



Large, blue
Art.No. 120050



Large, red
Art.No. 120055

Schine Pill Box is a storage box for medicins. The box is classified as a Med-tech product Class I.

- Practical storage of pills
- Easy opening and closing with transparent roller shutter system
- No unattached parts
- Large chambers with rounded corners
- Easy to transport
- Colored base to provide contrast toward the medicins.
- The large model also has prescription card in a plastic pouch on the back side (NB! The plastic pouch is not suitable for hand wash or dishwasher)

CARUMED

Företagshusvägen 5, 244 93 Kävlinge, Sweden
Tel: +46 (0) 46 - 73 41 00
info@carumed.se, www.carumed.se



Because we care!

Schine Pill Box

Small

Compartments: 8 (1-7 are labeled)
Weight: 60 grams
Length: 10,3 cms
Width: 5,4 cms
Height: 3,2 cms
Material: Recycable plastic
Washing: Can be hand washed

Large

Compartments: 28 (marked with numbers,
days and symbols)
Weight: 175 grams
Length: 19,7 cms
Width: 11,8 cms
Height: 3,1 cms
Material: Recycable plastic
Washing: Dishwasher top rack
(not the plastic pouch)




Företagshusvägen 5, 244 93 Kävlinge, Sweden
Tel: +46 (0) 46 - 73 41 00
info@carumed.se, www.carumed.se

CARUMED

Because we care!

Schine Pill Box



CARUMED

Företagshusvägen 5, 244 93 Kävlinge, Sweden

Tel: +46 (0) 46 - 73 41 00

info@carumed.se, www.carumed.se